



Instructions For Use

Easy One-Handed Operation Occludes Blood Flow In 5 Quick Steps

The **C-A-T**TM is Delivered in its One-Handed Configuration

This is the recommended carrying configuration



1

Insert the Wounded Extremity

Through the loop of the Self-Adhering Band



2

Pull the Self-Adhering Band Tight

And securely fasten the band back on itself



3

Adhere the Band Around the Arm

Do not adhere the band past the Windlass ClipTM



4

Twist the Windlass RodTM

Until bright red bleeding has stopped



5

Lock the Rod With the Windlass ClipTM

Bleeding is controlled



6

Adhere the Band Over the Rod

For small extremities, continue adhering the band around the extremity and over the Windlass RodTM



7

Secure the Rod and Band With the Strap

Adhering it to the opposite hook on the Windlass ClipTM.

The C-A-TTM is ready for transport



For added security and **always before moving a patient**, secure the Windlass RodTM with the Windlass StrapTM. For small extremities, also secure the Self-Adhering Band under the Windlass StrapTM.

Storing in the One-Handed Configuration

1

Pass 6" of the Self-Adhering Band

through the inside slit in the buckle, fold it back, and adhere the band to itself



2

Flatten the Loop of the C-A-TTM

Formed by the Self-Adhering Band. Placing the buckle in the middle



3

Fold the C-A-TTM in Half

Placing the buckle at one end





Instructions For Use

Application To A Leg Or Trapped Extremity



Combat Application Tourniquet™

1 Route Band Around the Leg

And pass the free-running end through the inside slit of the buckle



2 Pass Band Through the Outside Slit

of the buckle. This Utilizes the Friction Adaptor Buckle which will lock the band in place



3 Pull the Band Tight

And securely fasten the band back on itself



4 Twist the Windlass Rod™

Until bright red bleeding has stopped



5 Lock the Rod With the Clip



Bleeding is controlled

6 Secure the Rod With the Strap



C-A-T™ is ready for transport

Utilization of the Friction Adaptor Buckle

When using two hands to apply the C-A-T or if the Self-Adhering Band becomes exceptionally dirty, use the Friction Adaptor Buckle to lock the band in place

1 Pass Band Through the Inside Slit

Of the buckle



2 Pass Band Through the Outside Slit

Of the buckle. This will lock the band in place



3 Pull the Band Tight

And securely fasten the band back on itself



North American Rescue Products is the Exclusive U.S. Distributor of the C-A-T™.

www.NARescue.com

info@NARescue.com

888-689-6277

WARNING:

This product is intended for qualified medical professionals trained in the application of a tourniquet. It is not intended for the general public. Use only as directed by your EMS authority or under the supervision of a physician. Read entire Instructions For Use prior to using this product. Inappropriate application can result in serious injury or loss of limb. The use of any tourniquet is **A LAST RESORT** and should only be employed when bleeding cannot be stopped and the situation is life threatening.

LIMITATION OF LIABILITY:

Phil Durango, LLC, its employees, agents, contractors, suppliers, and distributors shall assume no liability for injury or damages arising from the application and use of the Combat Application Tourniquet™ (C-A-T™). The user assumes all risk of liability.

Phil Durango, LLC P 303.526.7847 F 303.736.4074 C-A-T@earthlink.net www.PhilDurango.com